



# Ambr<sup>®</sup> 250 Modular

Consumables and Accessories

Simplifying Progress

SARTORIUS

# Consumables and Accessories for the Ambr<sup>®</sup> 250 Modular System

Ambr<sup>®</sup> 250 Modular consumables are available to order individually or as part of a supply contract. Please contact your representative for more details.





# Vessels, Consumables and Accessories

## Ambr® 250 Modular Vessels

Bioreactor vessels designed with dual impellers either pitch-blade for mammalian or Rushton type for microbial cultures, configured as standard with integrated DO spot and pH probe. All bioreactors are irradiated and individually wrapped.

### Ambr® 250 “Easy-Connect” disposable bioreactors

| Part number | Description                                                             | Units per box |
|-------------|-------------------------------------------------------------------------|---------------|
| 001-2A23    | Mammalian vessel with dual pitch-blade impellers and baffles            | 2             |
| 001-2A13    | Microbial vessel with dual Rushton type impellers and baffles           | 2             |
| 001-2A33    | Mammalian low-shear vessel, w/o baffles, single ‘elephant ear’ impeller | 2             |

### Ambr® 250 Other Consumables and Accessories

| Part number | Description                     | Units per box |
|-------------|---------------------------------|---------------|
| 001-2A41    | Dust Cover                      | 1             |
| 001-2A42    | Light Protection Cover          | 1             |
| 001-2A52    | Reservoir Cleaning Kit          | 1             |
| 001-2A44    | Bioreactor Vessel Gas Filters   | 48            |
| 17575----K  | Minisart® SRP25 Syringe Filters | 50            |
| 17877---UPN | Midisart® 200 Filters           | 100           |
| 001-2A40    | Septum Caps                     | 8             |

### Ambr® 250 Modular Mammalian Vessel

Part no. 001-2A23

Ambr® 250 Modular bioreactor vessels, with baffles, dual pitch-blade impellers. Pack contains 2 irradiated bioreactors with caps, individually wrapped.



### Ambr® 250 Modular Microbial Vessel

Part no. 001-2A13

Ambr® 250 Modular bioreactor vessels, with baffles, dual Rushton type impellers. Pack contains 2 irradiated bioreactors with caps, individually wrapped.



### Ambr® 250 Modular Low Shear Vessel

Part no. 001-2A33

Ambr® 250 Modular bioreactor vessels, no baffles, single 'elephant ear' impeller. Pack contains 2 irradiated bioreactors with caps, individually wrapped.



# Covers

## Dust Cover

Part no. 001-2A41

Dust cover for Ambr® 250 Modular prevents dust from settling on the system.



## Light Protection Cover

Part no. 001-2A42

The light protection cover for Ambr® 250 Modular protects light sensitive substances inside the bioreactor vessel and its reservoirs from light.



# Gas Filters

## Ambr® 250 Modular Vessel Gas Filters

Part no. 001-2A44

Replacement gas filters for all Ambr® 250 Modular vessel variants. In case of a foam-out event, you can bring the bioreactor vessel under a biosafety cabinet and exchange the gas filters to save your experiment.

Pack contains 32 + 16 irradiated gas filters.



# Vessel Septum Caps

## Septum Caps Kit

Part no. 001-2A40

Every Ambr® 250 Modular vessel is equipped with a septum cap. In case you may want to exchange a septum cap on a vessel, you can use an irradiated septum cap of this kit.

Pack contains 8 irradiated septum caps, individually wrapped.




# Cleaning Kit

## Ambr® 250 Modular Cleaning Kit

Part no. 001-2A52

With the cleaning kit you have the right tools at hand to clean the reservoir drain of your Ambr® 250 Modular system. It contains an assortment of tools.



## Minisart® SRP25 Syringe Filters

Part no. 17575----K

Filters for system gas inlets. Unsterile, protects aeration module from particles.

Pack contains 50 unsterile gas filters.



## Midisart® 2000 Filters

Part no. 17877----UPN

Filters for system gas outlet. Unsterile, protects aeration module from particles.

Pack contains 100 unsterile gas filters.



## Germany

Sartorius Lab Instruments GmbH & Co. KG  
Otto-Brenner-Strasse 20  
37079 Goettingen  
Phone +49 551 308 0

## USA

Sartorius Corporation  
565 Johnson Avenue  
Bohemia, NY 11716  
Phone +1 631 254 4249  
Toll-free +1 800 635 2906

 For further contacts, visit  
[www.sartorius.com](http://www.sartorius.com)